

Putting Into Service

- 1 Switch off the TANGO controller.


- 2 Connect the 15-pin plug (1) and the HDI port of the TANGO controller using the connecting cable.
- 3 **Optional:** Connect an additional operating device to the 9-pin socket (2) (e. g. external Z-wheel-key).
- 4 Switch on the TANGO controller.
- 5 **Optional:** Configure the ERGODRIVE by means of SwitchBoard or TANGO instruction set.

Adjusting the Coaxial Drive

Vertical Position


Slide the drive knobs along the drive axis into the desired positions.

Horizontal Position


- 1 Loosen the 2 screws (1) using the Allen key.
- 2 Shift the coaxial drive (2) into the desired position.
- 3 Tighten the 2 screws (1) using the Allen key.

Adjusting the Drive Knobs' Ease of Movement

X-Drive Knob

- 1 Pull the X-drive knob (1) downwards.
- 2 Hold the X-drive knob (1) firmly.
- 3 Turn the adjustment ring (2) **clockwise to increase** or **counterclockwise to decrease** the resistance of the X-drive knob (1).

Y-Drive Knob

- 1 Hold the Y-drive knob (1) firmly.
- 2 Turn the adjustment ring (2) **clockwise to decrease** or **counterclockwise to increase** the resistance of the Y-drive knob (1).

Cleaning

WARNING

Danger to life from spreading hazardous substances.

In case of possible contamination with corrosive, infectious, toxic, radioactive or other harmful substances disinfect/decontaminate the ERGODRIVE professionally.

NOTICE

Risk of equipment damage from improper cleaning.

Ensure that no liquids enter the ERGODRIVE and especially the ports of the ERGODRIVE.

- 1 Switch off the TANGO controller.
- 2 Disconnect the ERGODRIVE from the TANGO controller.

- 3 Clean the ERGODRIVE with a lint-free cleaning cloth. If necessary, slightly dampen the cleaning cloth with water or a mild solvent-free disinfectant/cleaning agent.
- 4 Allow the ERGODRIVE to dry.
- 5 Put the ERGODRIVE into service again.

Maintenance

The ERGODRIVE is maintenance-free.

Taking Out of Service

- 1 Switch off the TANGO controller.
- 2 Disconnect the ERGODRIVE from the TANGO controller.

Service

You can reach our Customer Service Center Monday to Friday from 08:00 am to 04:00 pm CET. Please hold ready the information printed on the identification plate.

Tel.: +49 6441 9116-36

E-mail: service@marzhauser.com

Disposal

Referring to EU directive 2012/19/EU (WEEE) the ERGODRIVE is classified as electronic equipment and may not be disposed of with household waste. Märzhäuser Wetzlar will properly dispose of the ERGODRIVE. Also note the regulations outside the EU.

- 1 Take the ERGODRIVE out of service.
- 2 Send the ERGODRIVE to:

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Dimensions approx. 139 × 170 × 188 mm

Material aluminium, brass


Surface anodic coating, black lacquered

Declaration of Conformity

The ERGODRIVE was designed and produced in compliance with applicable directives and harmonized standards.

The Declaration of Conformity is available on request at: service@marzhauser.com

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Operating Manual ERGODRIVE

Operating Device for Motorized Microscope and Measuring Stages

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General Safety Instructions

The ERGODRIVE in combination with a Märzhäuser TANGO controller is intended and designed for manual movement of positioning systems.

- Read and understand this operating manual as well as the operating manual of the TANGO controller.
- Protect the ERGODRIVE from dust and liquids.
- The ERGODRIVE is not protected against corrosive, infectious, toxic, radioactive or other substances, which are harmful to health. When handling such substances, ensure that all legal requirements are fulfilled, especially the national regulations on accident prevention.
- Ensure that the ERGODRIVE, in conjunction with your application, meets the applicable safety regulations and legal requirements.
- The ERGODRIVE may only be repaired by our Customer Service Center.
- When lending or selling the device, make sure to provide this operating manual.

Delivery Contents


ERGODRIVE, connecting cable, Allen key (attached underneath the ERGODRIVE), operating manual

Configuration Options

The ERGODRIVE is individually configurable by means of SwitchBoard and TANGO instruction set:

- Configure travel speeds individually
- Move X-/Y-axes synchronously via X-drive knob
- Configure assignment of Z-wheel (any axis, LED 100)
- Enable/disable toggle mode for XY-/Z-key
- Configure F1-key (Snapshot, LED 100)
- Read states of all function keys
- Enable/disable noise filter
- Enable/disable ERGODRIVE

Function Description

No.	Description	Functions ¹
1	X-drive knob	move X-axis
2	Y-drive knob	move Y-axis
3	Z-wheel ²	move Z-axis

CAUTION

Avoid collisions between microscope stage and lens.

4	F1-key ³	not assigned
5	Z-key ^{2,3}	switch Z travel speed value 1 (preset): 0.1 mm/rev value 2: 1.0 mm/rev
6	XY-key ³	switch XY travel speed value 1 (preset): 1.0 mm/rev value 2: 14.0 mm/rev

¹ factory settings, functions individually configurable, refer to SwitchBoard user manual or TANGO instruction set

² only with ERGODRIVE 3

³ located on both sides of the ERGODRIVE

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Märzhäuser Wetzlar GmbH & Co. KG