

▶ Motorised microscope stage EK 75 x 50 Pilot

Märzhäuser Wetzlar GmbH & Co. KG

Single Cable Solution:

Communication between stage and controller plus electricity supply, all in one single, flexible cable (small cable diameter)!

equipped with USB converter, electricity supply via external wide-range power supply with corresponding adaptors (USA, UK, Australia, Europe)

One-Handed Operation:

Simply push slides into object holder and retrieve using one hand only to load and unload stage quickly.

integrated, very compact TANGO high-performance controller

Operation via integrated coaxial drive (up to 3-fold extendable) with 2 speed ranges for coarse positioning or fine adjustment. Torque setting for XY is infinitely variable possible.

