



► **Rotating and Centering Stage Z 150** with Object Guide AK 3

Diameter: 150 mm
(without vernier readings and binding screw)

Distance mounting platform to stage plate, top side: 14,5 mm

Weight: 780 g

Remark:

This rotating stage has a 360° graduation and two vernier readings.

The Z 150 may be centered right on the optical axis of a microscope.
The built-in stop has no effect on the lateral adjustment.

The extremely flat object guide AK 3 may be aligned on the stage to hold any standard size sample glass.

Order Information

Rotating and Centering Stage Z 150
(with object guide AK 3)
Order No.: 00-52-150-0000

Rotating and Centering Stage Z 150
Order No.: 00-52-100-0000

Accessory

For accessories please refer to:
www.marzhauser.com

Dimensions:

