


## ► Wafer Holders / Vacuum Wafer Holders for Upright Microscope- and Measuring Stages

### Wafer holder for 4" wafers (rotating)

Order no.:	00-29-440-8000
Dimensions (L x W x T):	116 x 116 x 13,2 mm
Adapter:	metal plate 116 x 116 mm (L x W)
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 100 x 100, SCAN <sup>plus</sup> 100 x 100


#### Dimensions



### Wafer holder for 5" and 6" wafers (rotating)

Order no.:	00-29-640-8000
Dimensions (L x W x T):	168 x 168 x 16,2 mm
Adapter:	metal plate 168 x 168 mm (L x W)
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 150 x 150, MTM 150 x 150 MR, SCAN 150 x 150


#### Dimensions



### Wafer holder for 5", 6" and 8" wafers (rotating)

Order no.:	00-29-840-8000
Dimensions (L x W x T):	246 x 246 x 16 mm
Adapter:	metal plate 246 x 246 mm (L x W)
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 200 x 200, MTM 200 x 200 MR, SCAN 200 x 200


#### Dimensions



### Wafer holder for 6", 8" and 12" wafers (rotating)

Order no.:	00-29-340-8000
Dimensions (D x T):	300,4 x 16 mm
Adapter:	metal plate 246 x 246 mm (L x W)
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 200 x 200, SCAN 200 x 200


#### Dimensions



### Vacuum wafer holder for 5" and 6" wafer (rotating)

Order no.:	00-29-650-8000
Dimensions (L x W x T):	168 x 168 x 24,5 mm
Adapter:	metal plate 168 x 168 mm (L x W)
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 150 x 150, MTM 150 x 150 MR, SCAN 150 x 150


#### Dimensions



### Vacuum wafer holder for 6" and 8" wafer (rotating)

Order no.:	00-29-850-8000
Dimensions (L x W x T):	200,4 x 200,4 x 24,5 mm
Adapter:	metal plate 246 x 246 mm (L x W)
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 200 x 200, SCAN 200 x 200

#### Dimensions



### Vacuum wafer holder for 8" and 12" wafer (rotating)

Order no.:	00-29-350-8000
Dimensions (D x T):	300,4 x 18,5 mm
Adapter:	3 threads in top plate of scanning stage
Material:	aluminium, steel and plastics
Surface:	anodic coating, black lacquered
Compatibility:	EK 300 x 300, SCAN 300 x 300

#### Dimensions

