


► Stage Inserts for inverse Microscope and Measuring Stages

Stage insert for 1 slide 3" x 1" (76 x 26 mm)

Order no.: 00-29-411-8010
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for 1 slide 3" x 1" (76 x 26 mm), rotating

Order no.: 90-29-411-8090
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for 2 slides 3" x 1" (76 x 26 mm)

Order no.: 00-29-411-8180
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for 3 slides 3" x 1" (76 x 26 mm)

Order no.: 38-29-411-8110
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for 4 slides 3" x 1" (76 x 26 mm)

Order no.: 00-29-411-8220
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Universal stage insert for samples Ø 20 - 74 mm and slides 3" x 1" (76 x 26 mm) and 3" x 2" (76 x 52 mm)

Order no.: 00-29-411-8150
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered

Dimensions




► Stage Inserts for inverse Microscope and Measuring Stages

Stage insert for Petri dishes, Ø 36 mm

Order no.: 00-29-411-8050
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for Petri dishes, Ø 54 mm

Order no.: 00-29-411-8060
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for Petri dishes, Ø 65 mm

Order no.: 00-29-411-8070
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for Petri dishes, Ø 88,5 mm

Order no.: 00-29-411-8080
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for 3 samples, Ø 41 mm

Order no.: 48-24-555-1000
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered


Dimensions



Stage insert for multi well plates

Order no.: 00-29-411-8020
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered

Dimensions



Stage insert for Terasaki 60/72/120

Order no.: 00-29-411-8030
Dimensions (L x W x T): 160 x 110 x 8 mm
Material: aluminium, steel
Surface: anodic coating, black lacquered

Dimensions

