



## 2-Axis Joystick for TANGO Controllers

### Specifications

Number of axes	2 (XY)
Interface	15-pin HD D-Sub plug connection to TANGO HDI
Material/surface	control panel: membrane keyboard housing: ABS-plastic with haptic coating housing bottom: stainless steel
Dimensions	229 x 122 x 90 mm (L x W x H)
Weight	approx. 1.0 kg

### Order Information

2-axis joystick  
 Order No.: 00-76-200-0820

### Product Features

- ▶ sensitive, manual positioning
- ▶ various configuration options
- ▶ free and precise programming of travel speed
- ▶ control of the LED illumination  
LED 100 possible
- ▶ 4 freely programmable function keys  
(e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

### Dimensions

