



## 4-Axis Joystick with Multi-Function Wheel for TANGO Controllers

### Specifications

Number of axes	4 (XYZA)
Interface	15-pin HD D-Sub plug connection to TANGO HDI
Material/surface	control panel: membrane keyboard housing: ABS-plastic with haptic coating housing bottom: stainless steel
Dimensions	229 x 139 x 96 mm (L x W x H)
Weight	approx. 1.0 kg

### Order Information

4-axis joystick  
 with Multi-Function Wheel  
 Order No.: 00-76-400-0820

### Product Features

- ▶ sensitive, manual positioning
- ▶ side-mounted multi-function wheel allows ergonomic operation
- ▶ various configuration options
- ▶ assignment of the multi-function wheel to any axis
- ▶ free and precise programming of travel speed
- ▶ control of the LED illumination LED 100 with multi-function wheel possible
- ▶ 4 freely programmable function keys (e. g. for position capture, velocity select etc.)
- ▶ membrane keyboard with high durability
- ▶ soft touch and insensitive haptic coating on the side surfaces
- ▶ easy cleaning

### Dimensions

